Addressed to:
AQUAMANA S.r.l.
Via S. Francesco di Sales 56
Carmagnola (TO), 10022

Note: this is an internal translation.

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Prot. n. L033/18 LAB

Object: verification of the containment of formaldehyde emission after an accidental spill, by using the products Novok-s (spray) and Novok-pow

As per your request, we evaluated the emission in question, after simulating an accidental spill of formaldehyde, a substance that is largely used in hospitals to fix biopsy specimens. The need for a containment of emissions, even accidental ones, derives from the urgency of reducing the exposure of operators to such substance (officially reclassified as carcinogenic - category IB - by IARC) as much as possible.

Considering these assumptions, the basic procedure to evaluate this containment was the following:

- a spill of 5 cm³ of 4% formaldehyde on a non-porous surface was simulated;
- the whole spill was covered with the powder “Novok-pow” (total height: about 3-6 mm of powder);
- Novok-s was applied directly on the surface of the powder (14 sprays);

In order to evaluate a potential performance improvement, a second test on formaldehyde vapours was carried out after removing the solidified mass; the removal was performed with the help of an absorbent cloth (14 sprays again). Finally, the whole area was cleaned.

The object of a third test was to evaluate the efficacy of containment using Novok-s only; it was carried out by following the procedure of the previous tests, omitting the use of Novok-pow.

All the tests were carried out under a switched off laminar flow hood, with the sampler about 1,5 m above the ground and 40 – 50 cm above the surface where the simulated spillage occurred. The hood was semi-open. We believe this procedure best simulates operators' working conditions when handling formalin. The choice of keeping the hood switched off implies a measurement by excess but guarantees more stable results, without introducing the suction variable; moreover, the data can be reasonably extrapolated in different but possible spill conditions, exactly like what occurs with a clear laboratory work table.

All the operations concerning the first test were performed in sequence, to make the data collection more realistic. Before proceeding with the third test, the test environment was cleaned and ventilated, by switching the hood on up to its maximum speed, for 10 minutes; the three measurements were repeated twice, in two different days. Both the initial environmental background and the maximum formaldehyde level after the spillage were measured in static conditions. The sampling temperature taken into consideration was between 22° and 24° C (room temperature).

The calibration of the used sampler allowed to collect air flows at 0,4 l/min, for a total volume of 4 litres of sampling, using Supelco cartridges (code 21026-U) – silica gel covered in 2,4-dinitrophenylhydrazine. The execution of the test follows the EPA method 0100 rev0 1996 when it comes to the collection and the EPA 8315 rev 1 dec 1996 when it comes to the analytical part.

The following tab shows the obtained results, including the abatement percentages that concern the spotted free emission value.

DATE	TEST DESCRIPTION	Formaldehyde concentration mg/m ³	AVERAGE VALUE mg/m ³	ABATEMENT PERCENTAGE
31/05/2018	Measurement of formaldehyde vapours emission from a spill of 5 cm ³	0.898		
31/05/2018	environmental background	0.012	0.013	
6/05/18	environmental background	0.014		
TEST 1				
31/05/2018	Measurement of formaldehyde vapours after covering the whole spill with Novok-Pow and spraying Novok-s 14 times on the surface of the mass	0.06	0.062	93.1
06/05/2018	Measurement of formaldehyde vapours after covering the whole spill with Novok-Pow and spraying Novok-s 14 times on the surface of the mass	0.063		
TEST 2				
31/05/2018	Measurement of formaldehyde vapours after TEST 1: cleaning procedure and further use of Novok-s (14 sprays)	0.027	0.029	96.7
06/05/2018	Measurement of formaldehyde vapours after TEST 1: cleaning procedure and further use of Novok-s (14 sprays)	0.032		
TEST 3				
31/05/2018	Measurement of formaldehyde vapours after covering the spill with Novok-s (14 sprays)	0.149	0.152	83.1
06/05/2018	Measurement of formaldehyde vapours after covering the spill with Novok-s (14 sprays)	0.155		

The measured values can be compared with the regulatory references and the national and international guidelines, summarised as follows:

LIMIT VALUES CONCERNING FORMALDEHYDE IN WORKING ENVIRONMENTS			
BODY	Limit value TWA- 8 hours	Limit value STEL- 15 minutes	Limit value CEILING
S.C.O.E.L.	0.369 mg/m ³	0.738 mg/m ³	
A.C.G.I.H./A.I.D.I.I.	0.12 mg/m ³		0.37 mg/m ³
LIMIT VALUES CONCERNING FORMALDEHYDE IN LIVING ENVIRONMENTS			
Ministry of Health, Circular no. 57 (22.06.1983) - within the living environments where products such as plywood, chipboards and cork mixtures are used		0.123 mg/m ³	
OMS - WHO: LOAEL "Lowest Observed Adverse Effect level" (30 minutes) - suggested limit		0.1 mg/m ³	

S.C.O.E.L.: “Scientific Committee on Occupational Exposure Limits” (UE)

A.C.G.I.H.: “American Conference of Governmental Industrial Hygienists”

A.I.D.I.I.: “Associazione italiana degli igenisti industriali” (Italian Association of Industrial Hygienists)

Considering all these assumptions, we can affirm that:

- 10 minutes after an accidental spill of 4% formaldehyde occurs, the operator's exposure level has to be considered critical and beyond the occupational limits suggested or required.
- The use of Novok-s, in combination with Novok-pow, brings the exposure back to a tolerable level, and below the levels considered tolerable for a non-working environment, with an abatement of more than 90%.
- A further cleaning procedure and other 14 sprays of Novok-s bring the emission values back to the environmental background, with an abatement of more than 95% compared to the initial value.
- The use of Novok-s only is effective to ensure an abatement level of the emission (beyond 80%) high enough to remain below the occupational exposure levels (S.C.O.E.L.), but above the ones applicable to a living environment.

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The accreditation certifies the technical competence of the laboratory limited to the scope detailed in the attached Enclosure. The scope may vary in the time. The management system requirements in ISO/IEC 17025:2005 (Section 4) are written in a language relevant to dei Laboratori di Prova operations and meet the principles of ISO 9001:2008 and are aligned with its pertinent requirements.

The present certificate is valid only if associated to the annexed schedule, and can be suspended or withdrawn at any time in the event of non fulfilment as ascertained by ACCREDIA.



The in force status of the accreditation may be checked in the WEB site (www.accredia.it) or on direct request to appointed Department.

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EUROLAB S.r.l.