


# 2-Lumen Cystometry Catheters


**Stericom's range of 2-lumen cystometry catheters offer the practitioner a choice of:**

- Pebax or PVC tubing
- Catheter length - 480mm, 550mm or 650mm
- Ch6 or Ch8 diameter
- Translucent or opaque tubing
- Bifurcation connector: plain, barrel-connector or 'Cath-Fix', for easy attachment to the patient




## **UDC 3776:**

- Ch6 Pebax catheter, cm markings on tubing
- Length of 2-lumen section: 480mm
- Total length from catheter tip to leur = 2,000mm



## **UDC 3716 (Ch6) UDC 3718 (Ch8):**

- PVC catheter, cm markings on tubing
- Colour-coded green (Ch6), blue (Ch8)
- Length of 2-lumen section: 550mm
- Total length from catheter tip to luer = 2,000mm



## **UDC 3726 (Ch6) UDC3728 (Ch8):**

- PVC catheter, cm markings on opaque 2-lumen tubing
- Colour-coded translucent green (Ch6), blue (Ch8)
- Length of 2-lumen section: 700mm
- Total length from catheter tip to leur = 2,000mm

**Contact us for full details**

**Pentland Medical Ltd. 157 Broomhall Crescent, Edinburgh, EH12 7PJ**

**Phone: 0131 467 5764 Fax: 0131 476 1574**

**E-Mail: [pentland.medical@virgin.net](mailto:pentland.medical@virgin.net) Web: [www.pentlandmedical.co.uk](http://www.pentlandmedical.co.uk)**