


Rectal/Abdominal Pressure Catheters

Stericom's range of rectal/abdominal catheters all feature:

- 2,000mm-long PVC tubing, for direct connection between the pressure transducer and the patient - translucent, red-tinted for easy line ID during the procedure
- Winged luer connectors, for ease of transducer connection

The practitioner is offered a choice of:

- PVC flap-valve - ideal for paediatrics
- Shaped, green PVC balloon
- 'Bovingdon' catheter with patented, shaped silicone balloon and introducer
- Catheters with or without stiffened introducer for finger-free insertion




UDC 3420/3421 Flap-valve Catheters

- Ch4.5 translucent, red-tinted tubing: length 2,000mm
- Winged luer - colour-coded red
- PVC flap-valve, open ended: ideal for paediatrics
- Either plain (UDC 3420) or with stiffened introducer (UDC 3421)
- Supplied in a box of 25x catheters



UDC 3430/3431 Green PVC Balloon Catheters

- Ch4.5 translucent, red-tinted tubing: length 2,000 mm
- Winged luer - colour coded red
- Green PVC balloon - with slit to assist irrigation and air expulsion
- Either plain (UDC 3430) or with stiffened introducer (UDC 3431)
- Supplied in a box of 25x catheters



UDC 3410/3414 Bovingdon Balloon Catheters

- Ch6 translucent, red-tinted tubing: length 2,000 mm
- Winged luer - colour-coded red
- Shaped silicone balloon - with slit to assist irrigation and air expulsion
- 100mm (UDC 3410) or 140mm (UDC 3414) flexible but stiffened introducer
- Supplied in a box of 25x catheters

Contact us for full details

Pentland Medical Ltd. 157 Broomhall Crescent, Edinburgh, EH12 7PJ

Phone: 0131 467 5764 Fax: 0131 476 1574

E-Mail: pentland.medical@virgin.net Web: www.pentlandmedical.co.uk