

# Single-Lumen Cystometry Catheters


## **Stericom's single-lumen cystometry catheters feature:**

- Single-lumen Ch4.5 cystometry catheter, in PVC
- Single-lumen Ch9 filling catheter, in PVC
- The catheters may be purchased as separate items, or together in a twin pack



### **UDC 3115 Single-Lumen Ch4.5 Cystometry Catheter**

- Ch4.5 PVC catheter, for filling/voiding pressures
- 2,000mm long for direct patient-transducer connection
- Translucent blue tubing
- Winged luer connector for secure & easy transducer connection
- Supplied in a box of 50x catheters



### **UDC 3210 2-Lumen Ch9 Filling Catheter**

- Ch9 PVC bladder filling catheter
- 400mm long
- Clear, uncoloured tubing
- Winged luer connector for secure & easy connection to filling/pump tube set
- Supplied in a box of 25x catheters



### **UDC 3312 Filling/Pressure Catheter Comination Pack**

- Contains 1x UDC 3115 & 1x UDC 3210 in each peel-pouch
- Ideal for 2-catheter filling/voiding cystometry measurements
- 'Piggy-back' catheterisation technique may be used
- Twin-pack saves space, inventory & cost
- Supplied in a box of 25x catheter sets



*Contact us for full details*

**Pentland Medical Ltd. 157 Broomhall Crescent, Edinburgh, EH12 7PJ**

**Phone: 0131 467 5764 Fax: 0131 476 1574**

**E-Mail: [pentland.medical@virgin.net](mailto:pentland.medical@virgin.net) Web: [www.pentlandmedical.co.uk](http://www.pentlandmedical.co.uk)**