


Furine® Vet

Fecal Management System



Video Tutorial

Furine Vet – the gentle choice
Product portfolio and accessories



Furine[®] Vet Fecal Management System

Clinical Purpose

Furine Vet is a non-invasive fecal management system (FMS) designed to safely divert fecal matter into a collection bag to protect the perianal skin from fecal exposure, maintain hygiene, reduce risk of infection, and improve patient care. The system is applied to the perianal skin of animal patients experiencing diarrhea or fecal incontinence.



Key Benefits

Furine Vet is beneficial for both patients and veterinary staff, as the product:

- Protects the patients perianal skin from diarrhea-related breakdown
- Reduces the risk of patient post-treatment trauma
- Helps minimize cross-contamination
- Enhances patient comfort
- Allows patient mobilization
- Reduces time spent on hygiene-related tasks
- Improves clinical efficiency
- Minimizes odor



ICU User Case: Canine Patient with Uncontrolled Diarrhea

Patient: A 10-year-old Schnauzer was admitted to the veterinary ICU with uncontrolled Diarrhea.

Clinical Problem: The patient could not control defecation whereby liquid feces repeatedly contaminate the perianal skin and tail base with risk of developing moisture-associated dermatitis. The staff had to clean the patient every 30–60 minutes.

Application: The Furine Vet FMS was applied to the perianal skin. The device adhered around the anus and directed feces into the collection bag.

Clinical Benefits: Feces were diverted away from the skin, preventing sore skin, dermatitis and maceration. The continuous fecal output was contained and measurable. Nursing time spent on repeated cleaning and trips outside for the dog were reduced notably. The patient remained cleaner, drier, and more comfortable.

Product Overview

Furine Vet

(ref. V104, V106, V108)

Furine Vet is a non-invasive fecal management system featuring a **dual silicone adhesive system** consisting of:

1. A soft and elastic adhesive base
2. A protective adhesive barrier film

The system diverts liquid stool into a dedicated collection bag, protecting the skin while maintaining a hygienic care environment.

Key features

- 100% silicone-based adhesive technology
- External application to the perianal skin
- Soft, flexible adhesive material to enhance patient comfort
- Hanger with secure locking mechanism
- Available in multiple sizes to accommodate different breeds and animal sizes



Furine Vet Collect

(ref. V205)

Furine Vet Collect is a collection bag designed for use with the Furine Vet FMS.

Key features:

- Opaque material for discreet handling
- Integrated viewing window for monitoring and measuring contents
- Secure closing mechanism to prevent leakage from the bag after use
- High performance charcoal filter allowing air release while containing odor



Securin Vet Skin Barrier

(ref. V305)

Securin Vet Skin Barrier is part of the dual silicone adhesive system used with the Furine Vet FMS. It functions as both a silicone adhesive and protective skin barrier film.

Key features:

- Forms a protective barrier on the skin
- Enhances secure attachment of Furine Vet to the skin
- Skin-friendly silicone based adhesive formulation
- Blue tint for improved visual guidance during application and removal



Securin Adhesive Remover Wipe

(ref. A001)

Securin Adhesive Remover Wipe is a siloxane-based medical cleansing wipe that is designed to prepare the skin prior to application of an adhesive product by removing excess oils. It also functions as an adhesive remover by dissolving and gently removing adhesive residues from the skin.

Key features:

- Functions as skin preparation wipe by removing excess oils from the skin
- Also serves as an adhesive remover wipe
- Gentle and non-stinging formulation
- Free from alcohol, oil, and perfume



Securin Adhesive Remover Spray

(ref. A011)

Securin Adhesive Remover Spray is a siloxane-based, gentle spray designed to support the safe and painless removal of adhesive products from the skin.

Key features:

- Supports gentle and painless removal of adhesive products
- Siloxane-based and sting-free formulation
- Free from alcohol, oil, and perfume
- Directional nozzle for precise application and reduced product waste
- Mechanical pump system – no propellants



Product:	Furine Vet			Furine Vet	Securin Vet	Securin Adhesive Remover	
	Small	Medium	Large	Collect	Skin Barrier	Wipe	Spray
Ref. No.:	V104	V106	V108	V205	V305	A001	A011

Contact info:

For further information, please visit www.furine.com or contact us at info@furine.com or by phone +45 7060 4131. Orders can be placed via order@furine.com

